

# VSI Reader ActionBar Installation & Use Cases

---

O.Burri, R. Guiet

BIOP, EPFL

Feb 2015

# Summary

---

- Installation
- Opening an Image
- Extracting ROIs Manually
- Defining Multiple ROIs
  - Manually
  - Automatically
  - As a Grid
- Batch Extract ROIs

# Description

Fast tiling microscopes are becoming commonplace but often rely on closed-source and closed-format systems to access the data.


This was the case with the Olympus Virtual Slide Scanner microscope. We needed a way to open and process the large images (>100k pixels) using the software in Fiji.

The result is a simple ActionBar that simplifies the most common tasks required when opening and dealing with these large datasets: ROI management and batch-processing.

## Note:

You will need Jerome Mutterer's excellent [ActionBar Plugin](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start) if you're going to use the VSI Reader.

Link: [http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action\\_bar:start](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start)




# Install VSI Reader ActionBar

You can install the VSI Reader via our [Fiji Update Site](http://fiji.sc/List_of_update_sites) (PTBIOP).

Link: [http://fiji.sc/List\\_of\\_update\\_sites](http://fiji.sc/List_of_update_sites)

Alternatively you can download it from the [BIOP Website](http://biop.epfl.ch/TOOL_VSI_Reader.html).

Link: [http://biop.epfl.ch/TOOL\\_VSI\\_Reader.html](http://biop.epfl.ch/TOOL_VSI_Reader.html)



[Start Update](#)

## Note:

You will need Jerome Mutterer's excellent [ActionBar Plugin](#) if you're going to use the VSI Reader.

Link: [http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action\\_bar:start](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start)


# Install VSI Reader ActionBar

You can install the VSI Reader via our [Fiji Update Site](http://fiji.sc/List_of_update_sites) (PTBIOP).

Link: [http://fiji.sc/List\\_of\\_update\\_sites](http://fiji.sc/List_of_update_sites)

Alternatively you can download it from the [BIOP Website](http://biop.epfl.ch/TOOL_VSI_Reader.html).

Link: [http://biop.epfl.ch/TOOL\\_VSI\\_Reader.html](http://biop.epfl.ch/TOOL_VSI_Reader.html)



## Note:

You will need Jerome Mutterer's excellent [ActionBar Plugin](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start) if you're going to use the VSI Reader.

Link: [http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action\\_bar:start](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start)



# Install VSI Reader ActionBar

You can install the VSI Reader via our [Fiji Update Site](http://fiji.sc/List_of_update_sites) (PTBIOP).


Link: [http://fiji.sc/List\\_of\\_update\\_sites](http://fiji.sc/List_of_update_sites)

Alternatively you can download it from the [BIOP Website](http://biop.epfl.ch/TOOL_VSI_Reader.html).

Link: [http://biop.epfl.ch/TOOL\\_VSI\\_Reader.html](http://biop.epfl.ch/TOOL_VSI_Reader.html)



Manage update sites



Select 'PTBIOP' and Close

Start Update

## Note:

You will need Jerome Mutterer's excellent [ActionBar Plugin](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start) if you're going to use the VSI Reader.

Link: [http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action\\_bar:start](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start)




# Install VSI Reader ActionBar

You can install the VSI Reader via our [Fiji Update Site](http://fiji.sc/List_of_update_sites) (PTBIOP).

Link: [http://fiji.sc/List\\_of\\_update\\_sites](http://fiji.sc/List_of_update_sites)

Alternatively you can download it from the [BIOP Website](http://biop.epfl.ch/TOOL_VSI_Reader.html).


Link: [http://biop.epfl.ch/TOOL\\_VSI\\_Reader.html](http://biop.epfl.ch/TOOL_VSI_Reader.html)




## Note:

You will need Jerome Mutterer's excellent [ActionBar Plugin](#) if you're going to use the VSI Reader.

Link: [http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action\\_bar:start](http://imagejdocu.tudor.lu/doku.php?id=plugin:utilities:action_bar:start)





# Start VSI Reader








# Start VSI Reader



# Open a VSI File 1/4





File Prompt pops-up




And Open

# Open a VSI File 2/4



Select the serie(s)  
you would like  
to visualize



A new window pops-up




And OK

(x,y) values are the size of  
the entire image in pixels


# Open a VSI File 3/4



Reduced version (Thumbnail)  
appears for each series






And a parameters window that contains  
information on the VSI File



Please, do not close this window,  
the macro needs information to process the  
images




# Open a VSI File 4/4



Original  
Size



Thumbnail  
Size

# Extract ROIs 1/4




If you define a ROI using the Rectangular Selection tool...


# Extract ROIs 2/4



And click on "Extract Current ROI"



A window pops-up asking for calibration information



Select:  
Correct Objective  
Correct camera


Should not be ticked if  
different magnifications  
are used

And OK



NB: Informations will be stored within the "Parameters" window, so you don't have to specify them for each new ROI in the series.




# Extract ROIs 3/4



ImageJ loads the raw pixel data...



...and the full size region is extracted







The unit is now 'microns' and the image is RGB.

If you specified "Imaging Mode" as **Fluorescence** Then you would get a Composite Image





# Extract ROIs 3/4




The unit is now 'microns'  
and the image is RGB.

If you specified  
"Imaging Mode" as Fluorescence  
Then you would get a Composite Image


# Extract ROIs 4/4





If you define a new ROI and extract



A new image pops-up




# Manual Define Multiple ROIs 1/3



Draw a ROI



The ROI is now in the ROI Manager





Informations  
(objective,  
imaging mode)  
have been  
already  
specified

Name	Date modified	Type
_Image01_	9/2/15 15:01	File folder
ROI Sets	9/2/15 15:01	File folder
Image_01.vsi	16/4/12 10:51	VSI File

Each Time you add a ROI, the ROIset is automatically saved  
in the ROI Sets folder


# Manual Define Multiple ROIs 2/3



If you draw another ROI



The new ROI is appended





Informations  
(objective,  
imaging mode)  
have been  
already  
specified

Name	Date modified	Type
_Image01_	9/2/15 15:01	File folder
ROI Sets	9/2/15 15:01	File folder
Image_01.vsi	16/4/12 10:51	VSI File

Each Time you add a ROI, the ROIset is automatically saved  
in the ROI Sets folder

# Manual Define Multiple ROIs 3/3




If you draw another ROI



And click on "Add to ROI manager"

The new ROI is appended



Select 'Show All' to see all the ROIs

Because the ROI set is automatically saved, you can do this for:



- Any number of ROIs in a Series
- Any number of Series per VSI File
- Any number of VSI Files

The only condition is that all VSI files be in the same folder




# Automatic Define Multiple ROIs 1/2


SIMPLE MODE



Select the folder where the VSI files are located



Click on "Batch Define ROIs"






Select either all series or the series you want to extract in each VSI file, separated by a comma

Select to use the whole image for each series.

# Automatic Define Multiple ROIs 2/2

SIMPLE MODE






Select Global Calibration to not be prompted for each series.

A ROI Set is saved for each series

PC > Swap (D:) > Experiment > ROI Sets >


Name	Date modified
7 - Series #0 (Thumbnail) from Image01.zip	11/2/15 15:16
14 - Series #1 (Thumbnail) from Image01.zip	11/2/15 15:16

# GRID Define Multiple ROIs






Select the size of the  
Individual images

Click on 'Make Grid From ROI'






# Extracting Multiple ROIs 1/2



You are prompted for the folder containing the VSI files  
(And the ROI Sets folder)





You can save all images in the same folder or  
in one folder per VSI file.

You can also downsample the images to  
speed up extraction

Click on "Process ROIs for extraction"

# Extracting Multiple ROIs 2/2



Within the folder containing the original image

Name	Type
_Image01_	File folder
Image01.vsi	Visual Studio Community Content Installer File

A new folder 'Extracted' is created

Name	Type
_Image01_	File folder
Image01.vsi	Visual Studio Community Content Installer File
Extracted	File folder

To save extracted ROI



As a 'tif' file

Image01.vsi\_Series\_0\_ROI\_1\_Scaling\_2\_RGB.tif

The VSI reader saves Brightfield images as RGB and Fluorescence images as composite TIFFs

# Fluorescent Acquisition Example









Image is a 3 Channel 16-bit  
Composite and will be saved  
as such

# End

- More is in preparation
  - Custom macros for defining ROIs
- Contact : mail at [Olivier and Romain](#)

- Acknowledgements
  - Anthony Grognez, Unité de Thérapie Cellulaire, CHUV
    - Widefield Dataset
  - Dr. Lalit Kumar Dubey, Harris Lab, EPFL
    - Fluorescent Dataset

The Fluorescent images should strictly not be used for any other purpose from you/or by any one in your organisation with out a written approval from Dr Lalit Kumar Dubey